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
COBAS AmpliPrep™ Instrument Maintenance

PURPOSE

The purpose of the SOP is to provide instructions for the frequency and performance of maintenance of the COBAS AmpliPrep Instrument.

SCOPE

This procedure is performed by properly trained NAT Lab personnel.


MATERIALS

- Deionized water
- Lint-free laboratory cloth
- 70% Isopropanol or an acceptable PCR Laboratory cleaner
- 10% sodium hypochlorite solution
- COBAS AmpliPrep wash cube container and System Wash Reagent
- COBAS AmpliPrep waste reservoir
- Tubing stands
- Replacement parts as appropriate
- TaqScreen WNV Test COBAS AmpliPrep Maintenance Log

PROCEDURE

AMPLILINK ASSISTED PROCEDURE

Step	Action
1	Verify that the Power ON LED light on the right side of the COBAS AmpliPrep is green.
2	Select the " Service " icon to display daily, weekly, and periodic service information.
3	Select the desired tab. <ul style="list-style-type: none"> 3.1 Due – Displays all tasks whose counters have passed the defined Elapsed %. 3.2 ALL – Displays all the service required for the TaqMan instrument. 3.3 Begin/End of Day (BOD) – Displays the service that will be performed when the next BOD occurs. 3.4 Service Log – Displays a log of completed service actions.
4	Select the desired service action to highlight
5	Select the <Perform> button to activate the Service Wizard (a wizard dialog box is displayed).
6	Perform the displayed step according to the displayed instructions and select the <Next> button when the displayed step is complete.


cobas AmpliPrep Maintenance Log

Serial Number: _____

Month/Year: _____

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31		
DAILY TASK																																	
Clean/Disinfect Spills																																	
Check/Replace Wash Buffer Cubes																																	
Empty Waste Reservoir																																	
Clean Rack Area																																	
Clean Initialization Posts																																	
Clean Sample Recognition Transceiver																																	
Check/Clean Tube Handler																																	
Check/Clean SPU Gripper																																	
Prime System (Normal Prime)																																	
Check Reagent Tip for Droplets																																	
Check Tubing for Leaks																																	
Clean Sensors and Light Barriers																																	
WEEKLY TASK																																	
Clean AmpliPrep Surfaces																																	
Clean SK24 Racks																																	
Clean SPU Racks																																	
Clean Reagent Racks																																	
Clean/Rinse Waste Containers																																	
Perform Extended Prime																																	
Not Used																																	
Out of Service																																	
Reviewer Initials																																	
PERIODIC TASK					Tech Initials				Reviewer Initials				PERIODIC TASK				Tech Initials				Reviewer Initials												
Replace Gripper O-rings (60 days)													Replace Filter Fan (180days)																				
Replace Reagent Tips (10,000 tests)													Change HEPA Filter (365 days)																				
Replace Syringe Plunger Tips (10,000 tests or 180 days)													Replace Syringes (730 days)																				

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